

# TruSight™ Tumor

Deep coverage of 26 genes for low mutant allele detection within regions involved in solid tumors delivered on proven next-generation sequencing technology.

## Highlights

- High accuracy, low mutant allele detection**  
Highly accurate variant analysis at limit of detection below 5% mutant allele frequency across 175 exonic regions with 1000x minimum coverage of each region
- Optimized for formalin-fixed, paraffin-embedded (FFPE) tissue**  
Exceptional sample success rates with minimal DNA input for accurate base calling even in degraded FFPE samples
- Best coverage of variants involved in solid tumors**  
Coverage of complete exons for analysis of molecular heterogeneity in highly relevant content selected from CAP and NCCN guidelines, and late stage clinical trials<sup>1-3</sup>


## Introduction

Next-generation sequencing (NGS) has helped dramatically increase knowledge of the pathways involved in cancer. Current tumor molecular profiling methods involve multiple single-gene assays, each targeting specific point mutations within hotspot regions. This iterative technique increases costs, requires more DNA and a longer time to generate data. In addition, single-gene genotyping assays reveal only a portion of what is happening within the tumor and can lead to a "one size fits all" view of biology. The same genomic information gained through these tests can be obtained faster and more cost-effectively using NGS technology, while providing additional information on emerging biomarkers implicated in recent clinical studies.

## High Accuracy with Low Mutant Allele Detection

Deep sequencing can reveal somatic variation in tumor sub-populations at levels never before possible. TruSight Tumor achieves limits of detection below 5% mutant allele frequency across 175 amplicons (Figure 1A). At this low level of detection, there is a risk that fixation artifacts are either classified as mutations or create an allele bias. This method leverages proprietary biochemistry and bioinformatics techniques to differentiate true variation from artifacts, yielding unprecedented accuracy and specificity (Figure 1B).

Figure 1: TruSight Tumor Performance



- TruSight Tumor accurately detected mutant allele frequency (MAF) below 5% on Horizon Discovery FFPE cell lines with titrated MAF variants for use in evaluating assay sensitivity.
- Over 200 libraries were generated from DNA extracted from various FFPE tissue types and laboratory sources. Sequencing of these libraries produced analytical sensitivity and specificity of > 99.9%, and a library success rate of > 99%.

## Optimized for Formalin-Fixed, Paraffin-Embedded (FFPE) Tissue

Formalin-fixation, paraffin-embedding techniques have long been the standard to protect tissues for downstream analysis and ease of archiving. This process introduces challenges for molecular approaches as it often degrades DNA into small fragments and has the potential to damage the DNA base pairs themselves. This library preparation and sequencing method was designed specifically to unlock the depth of molecular heterogeneity found in precious tumor samples (Figure 2). Illumina's optimized extraction and sample pre-qualification method (QC) provides maximum gDNA yield with archival FFPE tissues. This QC assay recommends a dilution factor based on sample quality, minimizing the DNA input into the library preparation assay resulting in highly uniform libraries (Figure 3).

## The Broadest Coverage of Content Implicated in Solid Tumors

The genes involved with solid tumors were carefully selected for this panel to include 26 genes and 175 amplicons of relevant content from CAP and NCCN guidelines, relevant publications, and late-phase pharmaceutical clinical trials. The TruSight Tumor content set provides coverage of entire exons in regions where variation has been cataloged in the COSMIC database<sup>4</sup> in oncogenes, and coverage of all exons in tumor suppressor genes (Figure 4). This provides a more comprehensive view of somatic variation in solid tumors including lung, colon, melanoma, gastric and ovarian—enabling clinical researchers to look beyond point mutations within hot spots in single genes.

Novel bioinformatics tools provide high levels of specificity at exceptionally low levels of mutant allele detection across the 26 genes and 175 amplicons provided (Table 2). The variant calling algorithms and filtering and annotation tools enable a scientist to easily analyze and understand the variation present in a sample.

Combined with Illumina's proven NGS technology, the TruSight Tumor content set delivers unmatched SNP, insertion, and deletion detection across the broadest regions of variation to enable a better understanding of underlying biology leading to improved clinical studies (Table 1 and 2).

## Manufactured for the Specific Needs of the Clinical Research Laboratory


TruSight content sets provides the option of higher quality standards, with a binding supply agreement to enable better laboratory management and minimize laboratory validations.

**Table 1: Confirmed Variants**

Gene	Variant(s) Confirmed	Reference Sequence	Alternate Sequence	Variant type	Frequency
<i>BRAF</i>	V600K	AC	TT	Dinucleotide	76.1%
<i>BRAF</i>	V600E	A	T	SNP	2.9%
<i>BRAF</i>	V600M	C	T	SNP	3.25%
<i>EGFR</i>	E746_A750del	AGGAATTAAAGAGAAGC	A	Deletion	4.6%
<i>EGFR</i>	L747_S752del_insV	GAATTAAGAGAAGCAACATC	GT	In/Del	18.0%
<i>EGFR</i>	delE746_S752insV	GAATTAAGAGAAGCAACATC	GT	In/Del	60.4%
<i>EGFR</i>	delL747_T751insP	ATTAAGAGAAGCAA	AC	In/Del	37.5%
<i>EGFR</i>	L858R	T	G	SNP	2.2%
<i>KIT</i>	S501insAY	C	CTGCCTA	Insertion	32.5%
<i>KIT</i>	V559del	GGTT	G	Deletion	27.5%
<i>KRAS</i>	G12D	C	T	SNP	4.9%
<i>KRAS</i>	G13D	C	T	SNP	10.3%
<i>KRAS</i>	D173D	A	G	SNP	16.5%
<i>PDGFRA</i>	D842_D845del	AGACATCATGCAT	A	Deletion	36.0%
<i>PIK3CA</i>	E545K	G	A	SNP	4.2%

Variants of note are confirmed in previously characterized FFPE samples with TruSight Tumor. This library prep method and bioinformatics solution enables detection of relevant complex variation including dinucleotide changes, insertions, and deletions.

**Figure 2: TruSight Tumor Workflow**



The TruSight Tumor content set is a comprehensive assay for examining the most relevant cancer genes involved in solid tumors, even in the most challenging FFPE samples. Researchers can go from extracted genomic DNA to data in two days with the MiSeq® system.



