



genes (exons only) plus exonic hotspots of an additional 39 genes, providing nearly 100% coverage of all targeted regions (Table 2).

This optimized oligo pool provides uniform coverage of the target regions, enabling > 500x coverage for > 95% of amplicons at > 5,000x mean coverage. This translates into 8 samples per run on a MiSeq® system with v3 reagents, providing the required sensitivity and accuracy to call rare variants with confidence.

### Simple, Integrated Workflow

The TruSight Myeloid Sequencing Panel offers a fully integrated DNA-to-data solution, including a streamlined workflow and automated data analysis with specific variant calling (Figure 1). Starting with 50 ng DNA isolated from blood, bone marrow, or fine needle aspirates (FNA), libraries are generated with the highly multiplexed oligonucleotide probes. Sample-specific indexes are added to each library, enabling pooling of libraries before sequencing for higher throughput. Pooled libraries are loaded into the MiSeq system for automated sequencing and data analysis. Analyzed data can be imported into VariantStudio software for accurate variant annotation, classification, and reporting. The entire process is completed in just three days.


### Optimized Assay Chemistry

The TruSight Myeloid assay begins with hybridizing a highly multiplexed pool of oligo pairs, one upstream and one downstream of each region of interest (Figure 2). Each oligo contains unique, target-specific sequence as well as a universal adapter sequence that is used in a subsequent amplification reaction. A proprietary extension-ligation reaction extends across the region of interest, followed by ligation to unite the two probes and yield a library of new templates with common ends. The resulting extension-ligation templates are PCR-amplified, incorporating two unique, library-specific indexes. Final reaction products are converted to a single-stranded, adapter-ligated normalized library using a bead-based protocol. The sequence-ready library can be loaded into the MiSeq reagent cartridge ready for sequencing on the MiSeq system without additional processing.

**Table 2: Coverage Details**

Cumulative target region size	~141 kb
Number of target genes	54
Amplicon size	~250 bp
Number of amplicons	568
Recommended mean coverage	5,000x
Target minimum coverage	500x
Percent exons covered at 500x	95

**Figure 1: Simple, Integrated Workflow**



TruSight Myeloid offers integrated library preparation and sequencing, as well as automated data analysis, creating a streamlined workflow that goes from DNA to data in 3 days.


### Data Analysis

Sequence data generated from TruSight Myeloid libraries are analyzed using the on-instrument MiSeq Reporter software. After demultiplexing and FASTQ file generation, the software uses a custom banded Smith-Waterman aligner to align the reads against the human hg19 reference genome to create BAM files. The Somatic Variant Caller then performs variant analysis for the specified regions. The outputs are VCF or gVCR files, which are text files that contain SNPs and small indels.

### Data Interpretation

VCF and gVCR files generated by the MiSeq Reporter software can be imported directly into the Illumina VariantStudio software application. This powerful tool enables clinical researchers to identify and classify disease-relevant variants quickly, and then communicate significant findings in concise, actionable reports. The intuitive framework offers flexible filtering options, streamlined variant classification, rapid and rich annotation, and customizable reporting options.

Figure 2: Optimized Assay Chemistry



The TruSight Myeloid assay enables simple, streamlined hybridization and amplification of targeted regions.

## Summary

The TruSight Myeloid Sequencing Panel enables clinical research laboratories to access expert-defined content for investigating genomic features associated with hematological malignancies. The optimized probe set provides comprehensive coverage of the regions known to be frequently mutated in myeloid cancers and myeloproliferative disorders. Using this panel, 54 genes can be analyzed in a single assay, saving time, money, and resources over current single-gene assays.

## Ordering Information

Product	Catalog No.
TruSight Myeloid Sequencing Panel (96 samples)	FC-130-1010

## Learn More

To learn more about the TruSight Myeloid Sequencing Panel and Illumina NGS technology, visit [www.illumina.com/trusightmyeloid](http://www.illumina.com/trusightmyeloid).

## References

1. Facts & Statistics. The Leukemia & Lymphoma Society. ([www.lls.org/#/diseaseinformation/getinformationsupport/factsstatistics/](http://www.lls.org/#/diseaseinformation/getinformationsupport/factsstatistics/))

AAAGAATGATAACAGTAAACACACTTCTGTTAACCTTAAGATTACTTGGATCCACTGATTCAACGTACCGTAAACGAACGTATCAATTGAGACTAAATATTAACGTACCATTAAAGAGCTACCGTCTTCTGTTAACCTTAAGATTACTTGGATCCACTGATTCA  
AATCAACGTACCGTAAACGAACGTATCAATTAAAGATTACTTGGATCCACTGATTCAACGTACCGTAAACGAACGTATCAATTGAGACTAAATATTAACGTACCATTAAAGAGCTACCGTCTTCTGTTAACCTTAAGATTACTTGGATCCACTGATTCA  
AACGACGAAAGAATGATAACAGTAAACACACTTCTGTTAACCTTAAGATTACTTGGATCCACTGATTCAACGTACCGTAAACGAACGTATCAATTGAGACTAAATATTAACGTACCATTAAAGAGCTACCGTCTTCTGTTAACCTTAAGATTACTTGGATCCACTGATTCA  
TTAAGGTACCATTAAGAGCTACCGTCAACAGTAAACACACTTCTGTTAACCTTAAGATTACTTGGATCCACTGATTCAACGTACCGTAAACGAACGTATCAATTGAGACTAAATATTAACGTACCATTAAAGAGCTACCGTCTTCTGTTAACCTTAAGATTACTTGGATCCACTGATTCA  
AAAGAATGATAACAGTAAACACACTTCTGTTAACCTTAAGATTACTTGGATCCACTGATTCAACGTACCGTAAACGAACGTATCAATTGAGACTAAATATTAACGTACCATTAAAGAGCTACCGTCTTCTGTTAACCTTAAGATTACTTGGATCCACTGATTCA  
AAGATTACTTGGATCCACTGATTCAACGTAAAGATTACTTGGATCCACTGATTCAACGTACCGTAAACGAACGTATCAATTGAGACTAAATATTAACGTACCATTAAAGAGCTACCGTCTTCTGTTAACCTTAAGATTACTTGGATCCACTGATTCA  
AACGTATCAATTGAGACTAAATATTAACGTACCATTAAAGATTACTTGGATCCACTGATTCAACGTACCGTAAACGAACGTATCAATTGAGACTAAATATTAACGTACCATTAAAGAGCTACCGTCTTCTGTTAACCTTAAGATTACTTGGATCCACTGATTCA

