

PRODUCT CODE	PRODUCT DESCRIPTION	QTY
RTK015PLD	CMV ELITE MGB™ Kit	100
STD015PLD	CMV ELITE Standard	8
RTS020PLD	EBV ELITE MGB™ Kit	100
STD020PLD	EBV ELITE Standard	16
CTR020PLD	EBV ELITE Positive Control	12
RTS031PLD	HSV1 ELITE MGB™ Kit	100
STD031PLD	HSV1 ELITE Standard	16
CTR031PLD	HSV1 ELITE Positive Control	12
RTS032PLD	HSV2 ELITE MGB™ Kit	100
STD032PLD	HSV2 ELITE Standard	16
CTR032PLD	HSV2 ELITE Positive Control	12
RTS035PLD	VZV ELITE MGB™ Kit	100
STD035PLD	VZV ELITE Standard	16
CTR035PLD	VZV ELITE Positive Control	12
RTS036PLD	HHV6 ELITE MGB™ Kit	100
STD036PLD	HHV6 ELITE Standard	16
CTR036PLD	HHV6 ELITE Positive Control	12
RTS038PLD	HHV8 ELITE MGB™ Kit	100
STD038PLD	HHV8 ELITE Standard	16
CTR038PLD	HHV8 ELITE Positive Control	12
RTS070PLD	Parvovirus B19 ELITE MGB™ Kit	100
STD070PLD	Parvovirus B19 ELITE Standard	16
CTR070PLD	Parvovirus B19 ELITE Positive Control	12
RTS078PLD	ADENOVIRUS ELITE MGB™ Kit	100
STD078PLD	ADENOVIRUS ELITE Standard	16
CTR078PLD	ADENOVIRUS ELITE Positive Control	12
RTS175PLD	BKV ELITE MGB™ Kit	100
STD175PLD	BKV ELITE Standard	16
CTR175PLD	BKV ELITE Positive Control	12


AVAILABLE SOON:

- ASPERGILLUS ELITE MGB™ Kit
- CHLAMYDIA TRACHOMATIS ELITE MGB™ Kit
- MTB ELITE MGB™ Kit

Please visit www.elitechgroup.com/elitechmgblegalnotice for complete licensing and warranty information.

- Leading edge assays with proprietary technology
- Higher sensitivity
- Monoreagent format
- Universal cycling conditions
- Melt-curve analysis

Distributor Contact Information :

Product not available in the US
Please contact your sales representative for product availability in your country

WORLDWIDE OFFICES:

Belgium T: +32 9 2820 531
Brazil T: +55 27 3025 1415
France T: +33 4 90 17 54 50
Italy T: +39 02 48 40 35 42

Middle East & Africa T: +971 4 375 2744
Serbia T: +381 11 2467119
Switzerland T: +41 26 663 86 60
The Netherlands T: +31 485 511 675

UK T: +44 1268 493377
United States T: +1 609 216 7360
855-ELITECH (855-354-8324)

Headquarters T: +33 1 41 45 07 13
www.elitechgroup.com
info@elitechgroup.com



Innovation in Real-Time PCR

Leading edge assays with proprietary technology

ELITE MGB™ is a revolutionary advance in Real-Time PCR chemistry. All ELITE MGB™ assays have been designed with our proprietary Minor Groove Binder (MGB) protein, Superbases™ and Eclipse® Dark Quencher technologies.

Our patented technology features shorter, overlapping probes that efficiently and accurately detect target DNA sequences, while offering greater sensitivity and specificity.

THE MGB protein is a synthetic molecule that binds to the minor groove of double stranded DNA molecules. In Real-Time PCR applications, MGB increases the stability of double stranded DNA complexes, specifically, the hybridization between the probe and the amplified DNA target. The increased DNA-DNA hybrid stability allows the design of shorter detection probes with higher specificity. Furthermore, the specificity of ELITE MGB™




probes increases the ability to discriminate between a perfectly matched sequence and a mismatched target compared to analogous MGB-free, longer counterpart probes.

SUPERBASES™ are specially engineered, nitrogen-based nucleotides. The proprietary design of these nucleotides minimizes secondary and tertiary structure observed in A-T or G-C rich DNA target sequences. This feature allows ELITE MGB™ assays to target even the most challenging gene sequences without compromising accuracy, clinical specificity or sensitivity.

The **ECLIPSE® DARK QUENCHER** is a proprietary fluorophore and dye quencher chemistry resulting in low background signals. Its key benefit is to ensure that every ELITE MGB™ assay will have the highest sensitivity by minimizing background signal interference.

ELITE MGB™ PROBES FEATURE HIGHER SENSITIVITY, SPECIFICITY AND LOW FLUORESCENCE BACKGROUND.

Universal Cycling Conditions



All the ELITE MGB™ assays utilize the same temperature and amplification profile. Known commonly as Universal Cycling Conditions, this profile consists of a three-step protocol: a denaturing step, an annealing step and an extension step. Universal Cycling Conditions enable laboratories to create multiple targets runs that can be run simultaneously on the same thermocycler.

DENATURATION STEP: The coiled ELITE MGB™ probe float freely in solution and will not emit fluorescence due to the close proximity of the fluorophore and dye-quencher.

ANNEALING STEP: The probe anneals to the target DNA sequence. In this conformation the fluorophore is separated from the quencher and emits light at a specific wavelength.

EXTENSION STEP: During DNA polymerization, the probe is displaced from the target gene sequence and reverts to its coiled-conformational state without being hydrolyzed. In this state, the probe is free in solution and will not emit fluorescence due to the close proximity of the fluorophore and dye-quencher.

Monoreagent format

ELITE MGB™ MASTERMIX has been reformulated to allow for a larger sample volume, significantly improving assay sensitivity. ELITE MGB™ mastermix is provided as a single solution in a ready-to-use format and includes all the components required for Real-Time PCR. No pipetting, mixing or reagent set-up is necessary; simply aliquot mastermix into the reaction wells and add your extracted sample.

Melt-curve analysis

THE MELT-CURVE ANALYSIS is a technique that monitors the relative fluorescence from the ELITE MGB™ probe as the temperature is gradually increased. This process creates a sequence specific dissociation profile, or melt-curve fingerprint, that is unique for a given double-stranded DNA fragment. For laboratories to utilize this valuable tool, Real-Time PCR methods must use DNA probes that are intact and undegraded. During Real-Time PCR, ELITE MGB™ Probes remain intact (unlike older Real-Time methods in which probes are degraded).


At low temperature, ELITE MGB™ probe is hybridized to the DNA target generating maximum fluorescence. Increasing the temperature causes the ELITE MGB™ probe to dissociate or “melt” from the DNA target. As the ELITE MGB™ probe “melts” the fluorophore and quencher are brought back in proximity resulting in lower fluorescence. See upper graph at left.

Creating a second graphic (first derivative) will identify the specific Melt Temperature (T_M) corresponding to the inflection point in the first graphic.

Because melt-curve analysis can resolve subtle pairing differences between the probe and target sequences, it is a valuable tool to identify clinically relevant DNA sequence alterations. A mismatch pairing could result in an underestimation of the quantity of bacterial or viral target present in the clinical sample. By analyzing the melt-curve profile at the end of each amplification, the accuracy, specificity and integrity of the quantitative analysis can be assessed.

As such, melt-curve analysis gives clinicians critical data that could be clinically or therapeutically important or impact the course of care for their patients.

Raw Fluorescence



First Derivative

